



Taq Plus DNA Polymerase

Cat. no. EC0201

Storage: -20°C

Concentration: 2.5 U/μl

Product Size

Product Components	EC0201	EC0202
Taq Plus DNA Polymerase	250 U	500 U
10× Taq Plus Buffer	1.8 ml	1.8 ml

PRODUCTOS BIO-LOGICOS

<http://www.pb-l.com.ar>

Introduction

Taq Plus DNA Polymerase is a special blend of *Taq* polymerase and *Pfu* Polymerase. It possesses both 5'-3' polymerase and exonuclease activity. The advantages of Taq Plus DNA Polymerase are high productivity and fidelity. Comparing to regular *Taq* polymerase, Taq Plus can efficiently amplify large DNA fragments (20Kb for simple templates and 10Kb for complex templates). Besides, it possesses higher extension rate and amplification efficiency than *Pfu* DNA Polymerase.

Taq Plus DNA Polymerase generates PCR products with 3'-dA overhangs that can be directly used in TA-cloning. To obtain higher cloning efficiency, however, PCR products could be purified and added 3-d'A overhangs before TA cloning procedures.

Unit Definition

One unit of Taq Plus DNA Polymerase is defined as the amount that incorporates 10 nmol of dNTPs into acid-insoluble material within 30 min at 74°C with activated salmon sperm DNA as the template-primer.

Storage Buffer

20 mM Tris-HCl (pH 8.0), 0.1 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50% glycerol

10× Taq Plus Buffer

200 mM Tris-HCl (pH 8.4), 200 mM KCl, 100 mM (NH₄)₂SO₄, 15 mM MgCl₂ and other components

- There are two types of 10× Taq Plus buffer can be chosen: Mg²⁺ plus and Mg²⁺ free.
- Mg²⁺ free buffer is supplied with separate 25 mM MgCl₂ Solution.
- Unless specifically requested, Mg²⁺ plus buffer will be supplied as regular component.

Applications

Amplify DNA fragments from complex templates (e.g. Genome) with high fidelity, for applications such as gene cloning, Site-directed mutagenesis, SNP etc.

Example

Note: The following example only for reference, user must set up optimal reaction system according to different reaction conditions such as different templates or primers etc.

1. To 50 μl PCR reaction system: 1kb fragment of human genomic DNA was amplified (if use different reaction system, please proportionally increase or decrease the amount of reaction components referring to this system).

Template	< 1 μg
Primer 1(10 μM)	1 μl
Primer 2(10 μM)	1 μl
10× Taq Plus Buffer	5 μl
dNTP Mixture (2.5 mM)	4 μl
Taq Plus (2.5 U/μl)	0.5-1 μl
ddH ₂ O	up to 50 μl

2. PCR cycle set-up:

94°C 3 min	} 30 cycles
94°C 30 sec	
55°C 30 sec	
72°C 1 min	
72°C 5 min	

3. Result detection: Load 5 μl PCR products to agarose gel for PCR detecting.